

FlexSys

The medical platform that adapts to your needs



FlexSys

**Modular. Great Value.
Made in Germany.**



Latest Technology - Made in Germany

FlexSys - the medical platform that adapts to your needs

Concept:

- Initially, customers purchase the base unit with their preferred module(s). Afterwards, they can expand their FlexSys with additional modules.
- Selection between five different laser or light modules.
- Up to three different modules can be included in one device.

Key Advantages:

- Save space and money:
The platform can be expanded and adapted to individual requirements without the necessity of buying another base unit.
- State-of-the-art technology:
The specifications and results of each module are at least comparable to competitive stand-alone systems.

Laser and Light Modules

FlexSys - add up to five modules to the base unit

EPL® UVB Module
308 nm



EPL® lite UVB Module
308 nm



Infrared Laser Module
1550 nm



Base Unit



Yellow Laser Module
577 nm



KTP Green Laser Module
532 nm

FlexSys - Purchasing Concept

The modular concept makes the FlexSys system highly economical. After the initial purchase you save the cost of the base unit, when acquiring additional modules.

Initial Purchase:

Base Unit + Module I



Yellow
Laser
Module

You purchase the base unit and your preferred module. The price is comparable to an equivalent stand-alone solution.

First Expansion:

Add Module II



Infrared Laser Module
1550 nm

Additional modules are purchased without the base unit. Prices are much lower compared to equivalent stand-alone systems.

Further Expansions:

Add Module III etc.




EPL® UVB Module
308 nm

You still only need to purchase the modules. Thus, the price is significantly lower than a comparable stand-alone system.


Indications

Many medical needs can be addressed



Yellow Laser module 577 nm

Gold standard wavelength - suitable for all skin types including dark skin types



Indications

- Rosacea
- Teleangiectasia
- Couperose
- Spider naevi
- Hemangioma
- Nevus flammeus (port-wine stain)
- Red spider veins
- Removal of benign pigmented lesions
- Removal of skin changes

Yellow Laser 577 nm

Powerful 5 W laser with handpiece and/or scanner

Facts & Features:

- HOPSL (high-power optically pumped semiconductor laser) with 5 W output.
- Spot size: 0.7 and 1.0 mm
- Optional scanner for fast and homogenous treatments of larger areas (up to 1 cm²).
- The scanner can also be used for single shots.



Base Unit



Yellow Laser Module
577 nm




FlexSys 577 nm | Post acne erythema

Settings: 18 J/cm², 80 % density, single session

Courtesy Dr Swapnil Shah, India

Pre-Treatment




After 1st Session



FlexSys 577 nm | Erythromelanosis follicularis

Settings: 20 mJ/cm², 20 ms, 80 % coverage

Courtesy Dr Swapnil Shah, India



FlexSys 577 nm | Granuloma Pyogenicum

Settings: 30 J/cm², 1 pass all along the lesion, 3 sessions once in 3 weeks

Courtesy Dr Swapnil Shah, India

Pre-Treatment



After 3rd Session



FlexSys 577 nm | Dermatosi s Papulosa Nigra

Settings: 24 J/cm², Single session

Courtesy Dr Swapnil Shah, India

Pre-Treatment



After 1st Session



FlexSys versus QuadroStarPRO YELLOW and Quanta 585

Comparison

GME - FlexSys with
577 nm module



Asclepion -
QuadroStarPRO
YELLOW



Quanta 585



FlexSys versus QuadroStarPRO YELLOW and Quanta 585

Comparison technical specifications

	FlexSys	QuadroStarPRO Yellow	Quanta 585
Wavelength	577 nm	577 nm	Quanta markets 585 nm, but states 580 nm \pm 6 nm, we measured one device: 577 nm
Power	5 Watt	5 Watt	5 Watt
Spot size handpiece	1 mm	1 mm	0.5; 1.0; 1.5 and 3 mm
Treatment area scanner	10 mm x 10 mm	15 mm x 15 mm	18 mm x 18 mm
Spot size scanner	0.7 mm	1 mm	1 mm
Skin cooling (scanner)	Air (optional)	Contact (integrated)	Air (optional)

FlexSys versus QuadroStarPRO YELLOW and Quanta 585

Comparison features and usability

Advantages of FlexSys

- Light source of all 3 lasers is similar but FlexSys has a superior scanner
- No need to change between scanner and handpiece
 - FlexSys scanner has a single shot mode for precise treatments
 - Switching between scanning and single shot mode by pressing a button on the scanner
- Smaller spot size of scanner for shorter pulses
 - Sensitive skin requires short pulses
 - 0.7 mm allows half the pulse duration compared to 1.0 mm
- Air cooling is better than contact cooling
 - No pressure on vessels = blood is not squeezed out
 - No reflection of light by glass plate = 100 % of laser energy reaching the skin

FlexSys versus QuadroStarPRO YELLOW and Quanta 585

Comparison of scanners

GME - FlexSys with
577 nm module



Asclepion -
QuadroStarPRO
YELLOW




Quanta 585



Cooling for FlexSys

Combination of FlexSys and Zimmer Cryo 7 / Cryo mini



1550 nm Laser for Wrinkles, Scars & Striae

Non-ablative fractional treatment: collagen formation without significant downtime



Indications

- Skin Resurfacing
- Lentigines (age spots)
- Lentigines solares
- Actinic keratosis
- Melasma
- Periorbital wrinkles
- Wrinkles
- Acne scars
- Surgical scars

Fractional 1550 nm Laser

Non-ablative fractional laser with scanner



Infrared Laser Module
1550 nm



Base Unit

Facts & Features:

- Power: up to 100 mJ per shot in less than 10 ms
- A state-of-the-art scanner enables you to choose the optimal shape, size and density for each treatment.

FlexSys 1550 nm | Dermatosi Papulosa Nigra + Scars

Settings: 80 / 70 / 60 mJ, 10 % / 8 % / 10 % Coverage, 3 sessions

Courtesy Dr Swapnil Shah, India

Pre-Treatment



After 3rd Session



FlexSys 1550 nm | Acne Scars

Settings: 70 mJ, 10 % Coverage, 3 sessions


Courtesy Dr Swapnil Shah, India



FlexSys 1550 nm | Erythematous Scars

Settings: 70 mJ, 8 % Coverage, 3 sessions

Courtesy Dr Swapnil Shah, India



FlexSys 1550 nm | Acne Scars

Settings: 70 mJ, 10 % Coverage, 3 Passes, 3 Sessions

Courtesy Dr Swapnil Shah, India



FlexSys 1550 nm | Acne and Scars

Settings: 60 mJ, 10 % Coverage, 3 Passes, 3 Sessions

Courtesy Dr Swapnil Shah, India

Pre-Treatment




After 3rd Session



FlexSys 1550 nm | Chicken Pox Scars

Settings: 1ST pass only on the lesions; 3 passes over the entire area, 2 Sessions

Courtesy Dr Swapnil Shah, India



FlexSys 1550 nm | Keloid

Settings: 3 Sessions by rotating the handpiece, 70 mJ, 10 % Coverage

Courtesy Dr Swapnil Shah, India

Pre-Treatment



After 2nd Session



FlexSys versus Icon / ResurFX / Re.pair

Comparison

GME - FlexSys with
1550 nm module



Cynosure - Icon



Lumenis - ResurFX



Solta Medical -
Fraxel Re:store



FlexSys versus Icon / ResurFX / Re.pair

Comparison technical specifications

	FlexSys	Icon	ResurFX	Re:store
Wavelength	1550 nm	1540 nm	1565 nm	1550 nm
Power	15 Watt	15 Watt	15 Watt	30 Watt
Technology	Scanner	Microlens	Scanner	Handpiece with roller tip
Spot size(s)	now: 0.35 mm future: 0.25 mm	0.25 mm and 0.13 mm	0.25 mm	0.14 mm
Max. energy per pulse	100 mJ	70 mJ	70 mJ	70 mJ
Pulse duration for 70 mJ	6 ms	15 ms	6 ms	3 ms
Cooling	Air cooling	No cooling	Contact cooling	Air cooling

FlexSys versus Icon / ResurFX / Re.pair

Comparison features and usability

Comparison to ResurFX by Lumenis:

- Very similar specs. Only difference is method of cooling
- ResurFX is much more expensive

Comparison to ICON by Cynosure:

- Icon uses micro-lenses: density is fixed and laser energy split into many individual pulses => density has to be achieved by overlapping and multiple passes
- ICON platform is very expensive

Comparison to Re:store by Solta:

- Re:store has smaller spots and shorter pulse durations. Unclear if this is of clinical advantage

EPL® UVB module 308 nm

Targeted phototherapy with monochromatic EPL® Excimer Pulsed Light



Indications

- Vitiligo
- Psoriasis
- Alopecia areata
- Atopic dermatitis
- Mycosis fungoides
- Lichen planus

EPL[®] - Excimer Pulsed Light

Monochromatic excimer lamp

EPL[®] UVB Module
308 nm



Base Unit

Facts & Features:


- Monochromatic 308 nm light generated by excimer lamp.
- Spot size: 17.5 cm²
- Adaptable to treatment areas of different sizes by interchangeable fused silica rods.
- Power density: up to 100 mW/cm²

FlexSys 308 nm | Vitiligo

Settings: Dose: 200 - 250 mJ/cm²


Courtesy Dr N Raboobee FFDerm (SA)

Pre-Treatment



N Raboobee

After 4th Session



N Raboobee

FlexSys 308 nm | Psoriasis

Settings: Dose: 600 - 900 mJ/cm²

Courtesy Dr N Raboobee FFDerm (SA)

Pre-Treatment



N Raboobee

After 4th Session



FlexSys 308 nm | Atopic Dermatitis


Settings: Dose: 300 - 400 mJ/cm²

Courtesy Dr N Raboobe FFDerm (SA)

Pre-Treatment



After 2nd Session



FlexSys 308 nm | Alopecia Areata

Settings: Dose: 200 - 600 mJ/cm²

Courtesy Dr N Raboobe FFDerm (SA)

Pre-Treatment



After 17th Session



FlexSys EPL versus Exciplex / 308 Excimer / KN-5000D

Comparison

GME - FlexSys with 308 nm module



Clarteis - Exciplex



Alma - 308 Excimer



Kernel - KN-5000D



FlexSys EPL versus Exciplex / 308 Excimer / KN-5000D

Comparison technical specifications

	FlexSys	Exciplex	308 Excimer	KN-5000D
Wavelength	308 nm	308 nm	308 nm	308 nm/311 nm
Power density	100 mW/cm ²	100 mW/cm ²	50 mW/cm ²	50 mW/cm ²
Max. treatment area	17.5 cm ²	25 cm ²	16 cm ²	20 cm ²
Time typical psoriasis dose (750 mJ/cm ²)	7 sec	7 sec	15 sec	15 sec
Time for typical vitiligo dose (300 mJ/cm ²)	3 sec	3 sec	6 sec	6 sec

FlexSys EPL versus Exciplex / 308 Excimer / KN-5000D

Comparison features and usability

Advantages of FlexSys

- Fused silica adapter for precise treatments with good visibility
- Large touch screen (of base unit) = customer friendly interface
- 100 mW/cm² for fast treatments (only Exciplex is similar)
- Handpiece is the most ergonomic one
- Option for intra-oral treatments
- Possibility to combine excimer lamp with laser modules in one platform

EPL[®] lite - Excimer Pulsed Light

EPL[®] lite 308 applicator: state-of-the art LED technology for occasional UVB treatments

EPL[®] lite UVB Module
308 nm




Base Unit

Facts & Features:

- Narrowband UVB light generated by LED.
- Spot size: 2.5 cm²
- Power density: up to 50 mW/cm²
- Fixed fused silica rod for precise treatment.
- Easy to use, economic, and very durable.

Green Laser: Vascular Lesions & Pigments

Wide range of applications with established technology - suitable for light skin types



Indications

- Rosacea
- Teleangiectasia
- Couperose
- Spider naevi
- Hemangioma
- Nevus flammeus (port-wine stain)
- Red spider veins
- Removal of benign pigmented lesions
- Removal of skin changes

“KTP” Green Laser 532 nm

Powerful 8 W laser with handpiece and/or scanner

Facts & Features:

- HOPSL (high-power optically pumped semiconductor laser) with 8 W output.
- Spot size: 1 mm
- Optional scanner for fast and homogenous treatments of larger areas (up to 1 cm²).
- The scanner can also be used for single shots.



Base Unit



KTP Green Laser Module
532 nm

Advantages



Unique Selling Points

Advantages of FlexSys

Versatility:	Wide range of treatment options for medical indications with 4 different wavelengths and 5 different modules.
Clinical effectiveness:	FlexSys modules are at least as effective as comparable stand-alone devices.
Cost efficiency:	The price of a FlexSys base unit with the first module is similar to a comparable stand-alone devices. Each additional modules is significantly less expensive.
Compact Design:	Small footprint for a wide range of treatment options.
Transportability:	The device can be used in multiple locations.
Made in Germany:	GME's high manufacturing standards guarantee excellent quality.



THANK YOU

Visit us at www.gmeonline.de